

	<p>Knee (PRO) Applicator: Cup Air flow: 2 Time: 1.0 - 3.0 min</p>		<p>Face (FACE) Applicator: Mask Air flow: 2 Time: 1.0 - 2.0 min</p>
	<p>Elbow (PRO) Applicator: Cup Air flow: 2 Time: 1.0 - 3.0 min</p>		<p>Forearm (PRO) Applicator: Cup Air flow: 2 Time: 1.0 - 3.0 min</p>
	<p>Shoulder (PRO) Applicator: Cup Air flow: 3 Time: 1.0 - 5.0 min</p>		<p>Skin disorder (SPOT) Applicator: Tube Air flow: 2 Time: 1.0 - 3.0 min</p>
	<p>Back (PRO) Applicator: Cup Air flow: 3 Time: 1.0 - 5.0 min</p>		<p>Jaw (PRO) Applicator: Cup Air flow: 2 Time: 1 min</p>
	<p>Ankle (PRO) Applicator: Cup Air flow: 2 Time: 1.0 - 3.0 min</p>		<p>Hip (PRO) Applicator: Cup Air flow: 3 Time: 1.0 - 5.0 min</p>
	<p>Wrist (PRO) Applicator: Cup/Tube Air flow: 2 Time: 1.0 - 2.0 min</p>		<p>Muscle pain (PRO) Applicator: Cup Air flow: 3 Time: 3.0 - 5.0 min</p>
	<p>Neck (PRO) Applicator: Cup Air flow: 2 Time: 1.0 - 3.0 min</p>		<p>Foot (PRO) Applicator: Cup Air flow: 2 Time: 1.0 - 2.0 min</p>

Note: all settings **can and should** be adjusted according the cold tolerance of the client as well as applicators **can be changed** depending on the area/problem.